



## PROTOCOL FOR IDENTIFICATION OF *C. JEJUNI*, *C. COLI* AND *C. LARI* BY GEL-BASED PCR

April 2021

Version 1

| History of changes |                    |  |            |
|--------------------|--------------------|--|------------|
| Version            | Sections changed   | Description of the change  | Date       |
| Version 1          |                    | Text edits. This protocol is aligned with the SOP for the interlaboratory study performed 2021 to validate the method. | 2021-04-21 |
| Draft              | New draft document | -  | 2020-11-02 |

# Content

- 1. Scope of the method .....2
- 2. References .....2
- 3. Protocol.....2
  - 3.1. DNA extraction .....2
  - 3.2. PCR set-up.....2
  - 3.3. Amplification.....3
  - 3.4. Electrophoresis .....4
  - 3.5. Results.....4

## 1. SCOPE OF THE METHOD

This protocol describes a gel-based multiplex PCR assay for confirmation and identification of *Campylobacter jejuni*, *C. coli* and *C. lari*. It also targets *C. upsaliensis*, but it is not as specific for this species. The method both detects the 23S rRNA of *Campylobacterales* and specific species targets.

## 2. REFERENCES

This protocol is based on ISO/CD 10272:2017/Amd1 [1], Wang et al. 2002 [2], and Chaban et al. 2009 [3].

## 3. PROTOCOL

### 3.1. DNA extraction

Transfer one sterile loop with approximately 1 µl colony material into 1 ml of 0.1xTE buffer (10 mM Tris and 0.1 mM EDTA, pH 8.0) and try to make a homogeneous suspension. Extract DNA with a thermal lysis step (15 min at 95 °C). After centrifugation for 3 min at 10 000 x g, the supernatant is used as DNA template. If the PCR analysis is not run on the same day, the template shall be stored at -20 °C.

### 3.2. PCR set-up

Prepare the master mix as described in Table 2 using the primers described in Table 1. The DNA template is used undiluted for PCR, but in case of only one band detected (species target or 23S target), dilute it 1:50 for repeated testing.

Table 1. Description of oligonucleotides and amplicons.

| Species (gene)                               | Primer | Sequence (5' – 3')                  | Amplicon size (bp) |
|--|--------|-------------------------------------|--------------------|
| <i>C. jejuni</i> ( <i>hipO</i> )             | CJF    | ACT TCT TTA TTG CTT GCT GC          | 323                |
|  | CJR    | GCC ACA ACA AGT AAA GAA GC          |                    |
| <i>C. coli</i> ( <i>glyA</i> )               | CCF    | GTA AAA CCA AAG CTT ATC GTG         | 126                |
|  | CCR    | TCC AGC AAT GTG TGC AAT G           |                    |
| <i>C. lari</i> ( <i>cpn60</i> ) <sup>a</sup> | JH0015 | TCT GCA AAT TCA GAT GAG AAA A       | 180                |
|  | JH0016 | TTT TTC AGT ATT TGT AAT GAA ATA TGG |                    |
| <i>C. upsaliensis</i> ( <i>glyA</i> )        | CUF    | AAT TGA AAC TCT TGC TAT CC          | 204                |
|  | CUR    | TCA TAC ATT TTA CCC GAG CT          |                    |
| <i>Campylobacterales</i> (23S rRNA)          | 23SF   | TAT ACC GGT AAG GAG TGC TGG AG      | 650                |
|  | 23SR   | ATC AAT TAA CCT TCG AGC ACC G       |                    |

<sup>a</sup> These primers, from reference [3], have replaced those in reference [2] since they detect both *C. lari* subsp. *lari* and *C. lari* subsp. *concheus*.

Table 2. Reagents.

| Reagent  | Final concentration   | Volume per sample (µl) |
|--|-----------------------|------------------------|
| Template DNA   | Maximum 250 ng        | 2,5                    |
| PCR grade water  | ---                   | As required            |
| PCR-buffer (without MgCl <sub>2</sub> ) <sup>a</sup>   | 1 x                   | As required            |
| MgCl <sub>2</sub> solution   | 2 mM                  | As required            |
| dNTP solution  | 0,2 mM of each dNTP   | As required            |
| PCR primers <i>C. jejuni</i> and <i>C. lari</i>  | 0,5 µM of each primer | As required            |
| PCR primers <i>C. coli</i>   | 1 µM of each primer   | As required            |
| PCR primers <i>C. upsaliensis</i>  | 2 µM of each primer   | As required            |
| PCR primers 23S rRNA   | 0,2 µM of each primer | As required            |
| <i>Taq</i> DNA polymerase  | 1,25 U                | As required            |
| Total volume   | ----                  | 25                     |
| <sup>a</sup> If the PCR buffer solution already contains MgCl <sub>2</sub> , the final concentration of MgCl <sub>2</sub> in the reaction mixture is adjusted to 2 mM. |                       |                        |

### 3.3. Amplification

Any well-maintained and calibrated cycler instrument can be used as long as it is appropriate for the method. Use the amplification program described in Table 3.

Table 3. Temperature-time program.

|   |             |
|---|-------------|
| Activation/initial denaturation <sup>a</sup>                              | 3 min/95 °C |
| Amplification   | 30 s/95 °C  |
|   | 30 s/59 °C  |
|   | 30 s/72 °C  |
| Number of cycles (amplification)  | 30          |
| Final extension   | 7 min/72 °C |
| <sup>a</sup> Use an initial denaturation time appropriate for the enzyme. |             |

### 3.4. Electrophoresis

The amplified PCR products are detected using a 1,5 % agarose gel. For example, the GeneRuler 100bp ready-to-use DNA ladder can be used for determination of size.

### 3.5. Results

The target sequences are detected if the sizes of the PCR product correspond to the expected length of the target DNA sequences (see Table 1). This should be determined using an appropriate DNA ladder and positive controls for each target. The detection of both species target and *Campylobacterales* 23S target are required for a positive result. If only one of the bands is detected, the analysis should be re-run with a 1:50 dilution of the template DNA. If still not detecting both targets, the sample is considered negative for *C. jejuni*, *C. coli*, *C. lari* and *C. upsaliensis*.

[1] ISO/CD 10272:2017/Amd1 “Amendment 1 of Microbiology of food and animal feeding stuffs – Horizontal method for detection and enumeration of *Campylobacter* spp. Part 1: detection method and Part 2: Colony-count technique”.

[2] Wang G, Clark CG, Taylor TM, Pucknell C, Barton C, Price L, Woodward DL, Rodgers FG. 2002. Colony multiplex PCR assay for identification and differentiation of *Campylobacter jejuni*, *C. coli*, *C. lari*, *C. upsaliensis*, and *C. fetus* subsp. *fetus*. *J Clin Microbiol.* 40(12):4744-4747

[3] Chaban B, Musil KM, Himsforth CG, Hill JE. 2009. Development of *cpn60*-based real-time quantitative PCR assays for the detection of 14 *Campylobacter* species and application to screening of canine fecal samples. *Appl Environ Microbiol.* 75(10):3055-61. doi: 10.1128/AEM.00101-09.